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
Tick length (sec)	4.5				
Full name	KNU(B)				
Online Name	KNU(B)	Test In Use	YES		
Test type	Photometric		LOW	HIGH	
		Test limit	*	*	Abs/min
Result unit	Abs/min	Initial absorbance	0.000	*	A
Number of Decim.	5	Dilution limit	*	*	Abs/min
		Secondary dil 1+	0.0	0.0	
		Critical limit	*	*	Abs/min
		Reflex test limit	*	*	Abs/min
		Reflex test			
Acceptance	Automatic	Reference class	LOW	HIGH	In Use
Dilution 1+	0.0				
Sample type	Sample type 5	Correction factor	1.00		
		Correction bias	0.00	Abs/min	
		Temperature	37.0	°C	
Calibration type	None				
Factor	1.00	Bias		0.00	
Bias correction in use	NO				
Manual QC in Use	NO	Routine QC in Use		NO	
Blank	None	Normal cuvette			
Reagent	KNU(B)-SUB	Volume(ul)		110	
Disp. With	Extra	Add. Volume(ul)		10	
Wash reagent	KNU(B)-SUB				
Reagent wash	Before dispense				
Syringe speed	Slow				
Incubation		Time (sec)		480	
Sample		Volume(ul)		30	
Disp. With	Extra	Add. Volume(ul)		50	
Dilution with	Water	Wash reagent		Water	
Incubation		Time (sec)		126	
Measurement	Kinetic				
Wavelength (nm)	405 nm	Side wavel. (nm)		None	
Curve type	Nonlinear				
Nonlinearity limits					
Conc. (Abs/min)	0.00000	Percent (%)		10	
Time (sec)	120	Points and intervals		7/18.0 (sec)	

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Reagent	KNU(B)-SUB	Lot	Expiry date (dd/mm/yyyy)
Stable on board (days)	1		
Alarm limit (ml)	1.0		
Information			
Vial volume	10 ml		
Barcode id			
Syringe speed	Slow		

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About Novozymes

Novozymes is the world leader in biological solutions. Together with customers, partners and the global community, we improve industrial performance while preserving the planet's resources and helping build better lives. As the world's largest provider of enzyme and microbial technologies, our bioinnovation enables higher agricultural yields, low-temperature washing, energy-efficient production, renewable fuel and many other benefits that we rely on today and in the future. We call it Rethink Tomorrow.

Luna no. 2020-09222-02

Novozymes A/S
Krogshøjvej 36
2880 Bagsvaerd
Denmark

Laws, regulations, and/or third-party rights may prevent customers from importing, using, processing, and/or reselling the products described herein in a given manner. Without separate, written agreement between the customer and Novozymes to such effect, information provided in this document "AS IS" is of good faith and does not constitute a representation or warranty of any kind and is subject to change without further notice.

novozymes.com