

Lipase Activity

PLU assay Determination by capillary-GC

C-GC analyses The injection volume is 0.7 μ l.
Carrier gas is helium, flow 1.7ml/min (46.2cm/sec), split 30(80°C).

Method Injector temperature 280°C.
Detector temperature 280°C.

Column temperature program:
0-2 min. 80°C.
2-5.4 min 80-250°C, temperature slope 50°C/min.
5.4-10 min 250°C.

After a series of injections the equipment is cleaned for 3 hours at the following temperatures:

Injector 350°C
Detector 350°C
Column 250°C.

For the C-GC analyses approx. 5 μ l of the reaction mixture is diluted with 995 μ l heptane. The C-GC yields all three components (1-propanol, propyllaurate and lauric acid) well separated. The retention times are approx. 0.9, 4.8, 6.3 minutes, respectively. All analyses are made in duplicate, and mixed standards are analyzed both before and after the samples. It is important, that all flasks, vials etc. are sealed immediately after use to prevent evaporation of propanol.

The following 6 mixed standards are used:

Table 1

Standard no.	Concentration (mM)		
	1-Propanol	Propyllaurate	Lauric acid
1	0	0	0
2	4	1	4
3	8	2	8
4	12	3	12
5	16	4	16
6	20	5	20



Calculations

From the analyses made on the standards a standard curve is made for each of the three materials. The slopes are used as response factors. The correlation coefficient is close to 1.

The conversion C of either lauric acid or 1-propanol into propyllaurate is calculated as:

$$C_{LA} = \frac{PL}{PL + LA} = \frac{\frac{\text{area PL}}{RF PL}}{\frac{\text{area PL}}{RF PL} + \frac{\text{area LA}}{RF LA}}$$

$$C_{PR} = \frac{PL}{PL + PR} = \frac{\frac{\text{area PL}}{RF PL}}{\frac{\text{area PL}}{RF PL} + \frac{\text{area PR}}{RF PR}}$$

Where

- C_{LA} : Fraction of lauric acid, which is transformed into propyllaurate.
 C_{PR} : Fraction of 1-propanol, which is transformed into propyllaurate.
PL: Concentration of propyllaurate (mM).
LA: Concentration of lauric acid (mM).
PR: Concentration of 1-propanol (mM).
Area PL: Area of propyllaurate peak on chromatogram (Counts).
RF PL: Response factor of propyllaurate from standard curve (Counts/mM).

The conversion is then used in the formula below to calculate the activity of the enzyme analyzed in PLU/g.

$$\text{Activity PLU/g} = \frac{M \times C}{W \times t}$$

Where

- M: Initial amount of 1-propanol and lauric acid, i.e. 40000 μ mole.
C: Fraction ester from the formula above.
W: Amount of dry matter of catalyst in g.
t: Reaction time in min.

Enzyme Business

Novo Nordisk A/S
Novo Allé
2880 Bagsvaerd
Denmark

Tel. +45 4444 8888
Fax +45 4444 1021
Telex 37560
enzymes@novo.dk
www.novo.dk/enzymes

Laws, regulations and third party rights may prevent customers from importing, processing, applying and/or reselling certain products in a given manner. It is the responsibility of the customers that their specific use of products from Novo Nordisk does not infringe relevant laws and regulations and, furthermore, does not infringe patents or other third party rights.

The contents of this document are subject to change without further notice.